

Product Datasheet

Human Interleukin Enhancer-binding Factor 3 ELISA Kit (orb1213099)

Description	Human Interleukin Enhancer-binding Factor 3, ILF3 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	ILF3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	95,338
Uniprot ID	Q12906
Solubility (25°C)	2400ng/L
Sensitivity	4.84ng/L
Entrez	3609
Expiration Date	Please enquire.