

Product Datasheet

Human Interferon-induced Gtp-binding Protein Mx2 ELISA Kit (orb1213100)

Description	Human Interferon-induced Gtp-binding Protein Mx2, MX2 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	MX2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	82,089
Uniprot ID	P20592
Solubility (25°C)	4800ng/L
Sensitivity	9.96ng/L
Entrez	4600
Expiration Date	Please enquire.