

Product Datasheet

Human Interferon-induced Gtp-binding Protein Mx1 ELISA Kit (orb1213101)

Description	Human Interferon-induced Gtp-binding Protein Mx1, MX1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	MX1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	75,520
Uniprot ID	P20591
Solubility (25°C)	2400ng/L
Sensitivity	5.23ng/L
Entrez	4599
Expiration Date	Please enquire.