

Product Datasheet

Human Interferon Epsilon ELISA Kit (orb1213103)

Description	Human Interferon Epsilon, IFNE ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	IFNE
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,414
Uniprot ID	Q86WN2
Solubility (25°C)	2560ng/L
Sensitivity	5.04ng/L
Entrez	338376
Expiration Date	Please enquire.