

Product Datasheet

Human Interferon Alpha-inducible Protein 6 ELISA Kit (orb1213104)

Description	Human Interferon Alpha-inducible Protein 6, IFI6 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	IFI6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	12,927
Uniprot ID	P09912
Solubility (25°C)	3200ng/L
Sensitivity	7.68ng/L
Entrez	2537
Expiration Date	Please enquire.