

Product Datasheet

Human Inter-alpha-trypsin Inhibitor Heavy Chain H2 ELISA Kit (orb1213107)

Description	Human Inter-alpha-trypsin Inhibitor Heavy Chain H2, ITIH2 ELISA Kit
Reactivity	Human
Range	1-400ug/mL
Tested Applications	ELISA
Target	ITIH2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	106,463
Uniprot ID	P19823
Solubility (25°C)	480µg/mL
Sensitivity	0.54µg/mL
Entrez	3698
Expiration Date	Please enquire.