

Product Datasheet

Human Integrin Beta-8 ELISA Kit (orb1213109)

Description	Human Integrin Beta-8, ITGB8 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	ITGB8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	85,632
Uniprot ID	P26012
Solubility (25°C)	48ng/ml
Sensitivity	0.053ng/ml
Entrez	3696
Expiration Date	Please enquire.