

Product Datasheet

Human Integrin Beta-7 ELISA Kit (orb1213110)

Description	Human Integrin Beta-7, ITGB7 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	ITGB7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	86,903
Uniprot ID	P26010
Solubility (25°C)	48ng/ml
Sensitivity	0.051ng/ml
Entrez	3695
Expiration Date	Please enquire.