

Product Datasheet

Human Insulin-like Growth Factor Binding Protein 5 ELISA Kit (orb1213118)

Description	Human Insulin-like Growth Factor Binding Protein 5, IGFBP-5 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	IGFBP5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,570
Uniprot ID	P24593
Solubility (25°C)	24ng/ml
Sensitivity	0.023ng/ml
Entrez	3488
Expiration Date	Please enquire.