

Product Datasheet

Human Insulin-like Growth Factor 2 MRNA-binding Protein 3 ELISA Kit (orb1213119)

Description	Human Insulin-like Growth Factor 2 MRNA-binding Protein 3, IGF2BP3 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	IGF2BP3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	63,705
Uniprot ID	000425
Solubility (25°C)	4800ng/L
Sensitivity	9.87ng/L
Entrez	10643
Expiration Date	Please enquire.