

Product Datasheet

Human Inositol 1,4 ELISA Kit (orb1213122)

Description	Human Inositol 1,4, 5-Trisphosphate Receptor Type 1, ITPR1 ELISA Kit
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	ITPR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	313,929
Uniprot ID	Q14643
Solubility (25°C)	960ng/L
Sensitivity	1.46ng/L
Entrez	3708
Expiration Date	Please enquire.