

Product Datasheet

Human Inhibitor of Growth Protein 4 ELISA Kit (orb1213123)

Description	Human Inhibitor of Growth Protein 4, ING4 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	ING4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,530
Uniprot ID	Q9UNL4
Solubility (25°C)	6400ng/L
Sensitivity	9.82ng/L
Entrez	51147
Expiration Date	Please enquire.