

Product Datasheet

Human Indoleamine 2 ELISA Kit (orb1213125)

Description	Human Indoleamine 2, 3-dioxygenase 2, IDO2 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	IDO2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,075
Uniprot ID	Q6ZQW0
Solubility (25°C)	1600ng/L
Sensitivity	3.62ng/L
Entrez	169355
Expiration Date	Please enquire.