

Product Datasheet

Human Inactive Tyrosine-protein Kinase 7 ELISA Kit (orb1213126)

Description	Human Inactive Tyrosine-protein Kinase 7, PTK7 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	CCK4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	118,392
Uniprot ID	Q13308
Solubility (25°C)	1600ng/L
Sensitivity	3.35ng/L
Entrez	5754
Expiration Date	Please enquire.