

Product Datasheet

Human Hyaluronidase-3 ELISA Kit (orb1213131)

Description	Human Hyaluronidase-3, HYAL3 ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	HYAL3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,501
Uniprot ID	043820
Solubility (25°C)	96ng/mL
Sensitivity	0.15ng/mL
Entrez	8372
Expiration Date	Please enquire.