

Product Datasheet

Human Hla Class 2 Histocompatibility Antigen ELISA Kit (orb1213138)

Description	Human Hla Class 2 Histocompatibility Antigen, Dr Alpha Chain, HLA-DRA ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	HLA-DRA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,621
Uniprot ID	P01903
Solubility (25°C)	2560ng/L
Sensitivity	5.48ng/L
Entrez	3122
Expiration Date	Please enquire.