

Product Datasheet

Human Histone Deacetylase 9 ELISA Kit (orb1213146)

Description	Human Histone Deacetylase 9, HDAC9 ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	HDAC9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	111,297
Uniprot ID	Q9UKV0
Solubility (25°C)	25.6ng/mL
Sensitivity	0.044ng/mL
Entrez	9734
Expiration Date	Please enquire.