

Product Datasheet

Human Histatin-3 ELISA Kit (orb1213148)

Description	Human Histatin-3, HTN3 ELISA Kit
Reactivity	Human
Range	0.05-20ug/mL
Tested Applications	ELISA
Target	HTN3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	6,149
Uniprot ID	P15516
Solubility (25°C)	24µg/mL
Sensitivity	0.029µg/mL
Entrez	3347
Expiration Date	Please enquire.