

Product Datasheet

Human Histatin-1 ELISA Kit (orb1213149)

Description	Human Histatin-1, HTN1 ELISA Kit
Reactivity	Human
Range	10-2000ng/mL
Tested Applications	ELISA
Target	HTN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	6,963
Uniprot ID	P15515
Solubility (25°C)	2400ng/mL
Sensitivity	4.39ng/mL
Entrez	3346
Expiration Date	Please enquire.