

Product Datasheet

Human Heat Shock Protein Hsp 90-beta ELISA Kit (orb1213160)

Description	Human Heat Shock Protein Hsp 90-beta, HSP90AB1 ELISA Kit
Reactivity	Human
Range	0.5-200ng/mL
Tested Applications	ELISA
Target	HSP90AB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	83,264
Uniprot ID	P08238
Solubility (25°C)	256ng/mL
Sensitivity	0.36ng/mL
Entrez	3326
Expiration Date	Please enquire.