

Product Datasheet

Human Gtp Cyclohydrolase 1 ELISA Kit (orb1213164)

Description	Human Gtp Cyclohydrolase 1, GCH1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	GCH1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,903
Uniprot ID	P30793
Solubility (25°C)	3200ng/L
Sensitivity	7.47ng/L
Entrez	2643
Expiration Date	Please enquire.