

Product Datasheet

Human G-protein Coupled Estrogen Receptor 1 ELISA Kit (orb1213170)

Description Human G-protein Coupled
Estrogen Receptor 1, GPR30
ELISA Kit

Reactivity Human

Range 0.5-200ng/L

**Tested
Applications** ELISA

Target GPER1

Storage 2-8°C

Note For research use only

Application notes Sandwich

Sample Types Serum, plasma, cell culture
supernates

MW 42,248

Uniprot ID [Q99527](#)

Solubility (25°C) 256ng/L

Sensitivity 0.26ng/L

Entrez [2852](#)

Expiration Date Please enquire.