

Product Datasheet

Human Glutathione Peroxidase 2 ELISA Kit (orb1213174)

Description	Human Glutathione Peroxidase 2, GPX2 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	GPX2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,954
Uniprot ID	P18283
Solubility (25°C)	3200ng/L
Sensitivity	7.73ng/L
Entrez	2877
Expiration Date	Please enquire.