

Product Datasheet

Human Glutamate Dehydrogenase 1 ELISA Kit (orb1213178)

Description	Human Glutamate Dehydrogenase 1, Mitochondrial, GLUD1 ELISA Kit
Reactivity	Human
Range	5-1000uU/mL
Tested Applications	ELISA
Target	GLUD1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	61,398
Uniprot ID	P00367
Solubility (25°C)	1280μU/mL
Sensitivity	2.38μU/mL
Entrez	2746
Expiration Date	Please enquire.