

Product Datasheet

Human Gastrotropin ELISA Kit (orb1213185)

Description	Human Gastrotropin, FABP6 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	FABP6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,371
Uniprot ID	P51161
Solubility (25°C)	48ng/ml
Sensitivity	0.048ng/ml
Entrez	2172
Expiration Date	Please enquire.