

## Product Datasheet

# Human G Protein-coupled Bile Acid Receptor 1 ELISA Kit (orb1213193)

|                            |                                                                |
|----------------------------|----------------------------------------------------------------|
| <b>Description</b>         | Human G Protein-coupled Bile Acid Receptor 1, GPBAR1 ELISA Kit |
| <b>Reactivity</b>          | Human                                                          |
| <b>Range</b>               | 0.05-20ng/mL                                                   |
| <b>Tested Applications</b> | ELISA                                                          |
| <b>Target</b>              | GPBAR1                                                         |
| <b>Storage</b>             | 2-8°C                                                          |
| <b>Note</b>                | For research use only                                          |
| <b>Application notes</b>   | Sandwich                                                       |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                         |
| <b>MW</b>                  | 35,248                                                         |
| <b>Uniprot ID</b>          | <b>Q8TDU6</b>                                                  |
| <b>Solubility (25°C)</b>   | 25.6ng/mL                                                      |
| <b>Sensitivity</b>         | 0.027ng/mL                                                     |
| <b>Entrez</b>              | <b>151306</b>                                                  |
| <b>Expiration Date</b>     | Please enquire.                                                |