

Product Datasheet

Human Fibroblast Growth Factor-13 ELISA Kit (orb1213200)

Description	Human Fibroblast Growth Factor-13, FGF-13 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	FGF13
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,564
Uniprot ID	Q92913
Solubility (25°C)	1600ng/L
Sensitivity	4.04ng/L
Entrez	2258
Expiration Date	Please enquire.