

Product Datasheet

Human Fibroblast Growth Factor Receptor 4 ELISA Kit (orb1213201)

Description	Human Fibroblast Growth Factor Receptor 4, FGFR4 ELISA Kit
Reactivity	Human
Range	1-400ng/mL
Tested Applications	ELISA
Target	FGFR4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	87,954
Uniprot ID	P22455
Solubility (25°C)	480ng/mL
Sensitivity	0.57ng/mL
Entrez	2264
Expiration Date	Please enquire.