

Product Datasheet

Human Fibroblast Growth Factor Receptor 3 ELISA Kit (orb1213202)

Description	Human Fibroblast Growth Factor Receptor 3, FGFR3 ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	FGFR3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	87,710
Uniprot ID	P22607
Solubility (25°C)	9600ng/L
Sensitivity	14.93ng/L
Entrez	2261
Expiration Date	Please enquire.