

Product Datasheet

Human Fibroblast Growth Factor 5 ELISA Kit (orb1213204)

Description	Human Fibroblast Growth Factor 5, FGF5 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	FGF5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,551
Uniprot ID	P12034
Solubility (25°C)	2560ng/L
Sensitivity	5.26ng/L
Entrez	2250
Expiration Date	Please enquire.