

Product Datasheet

Human Fibroblast Growth Factor 12 ELISA Kit (orb1213206)

Description	Human Fibroblast Growth Factor 12, FGF12 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	FGF12
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,399
Uniprot ID	P61328
Solubility (25°C)	4800ng/L
Sensitivity	10.31ng/L
Entrez	2257
Expiration Date	Please enquire.