

Product Datasheet

Human Fibroblast Growth Factor 10 ELISA Kit (orb1213207)

Description	Human Fibroblast Growth Factor 10, FGF10 ELISA Kit
Reactivity	Human
Range	1-400ng/L
Tested Applications	ELISA
Target	FGF10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,436
Uniprot ID	O15520
Solubility (25°C)	480ng/L
Sensitivity	0.57ng/L
Entrez	2255
Expiration Date	Please enquire.