

Product Datasheet

Human Fibrillin-3 ELISA Kit (orb1213208)

Description	Human Fibrillin-3, FBN3 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	MATN3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,817
Uniprot ID	O15232
Solubility (25°C)	48ng/ml
Sensitivity	0.046ng/ml
Entrez	4148
Expiration Date	Please enquire.