

Product Datasheet

Human Fatty Acid-binding Protein ELISA Kit (orb1213212)

Description	Human Fatty Acid-binding Protein, Brain, FABP7/BLBP/FABPB/MRG ELISA Kit
Reactivity	Human
Range	0.5-200ng/ml
Tested Applications	ELISA
Target	FABP7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,889
Uniprot ID	O15540
Solubility (25°C)	256ng/ml
Sensitivity	0.29ng/ml
Entrez	2173
Expiration Date	Please enquire.