

Product Datasheet

Human Extracellular Matrix Protein 1 ELISA Kit (orb1213216)

Description	Human Extracellular Matrix Protein 1, ECM1 ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	ECM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	60,674
Uniprot ID	Q16610
Solubility (25°C)	9600ng/L
Sensitivity	15.23ng/L
Entrez	1893
Expiration Date	Please enquire.