

Product Datasheet

Human Epithelial Discoidin Domain-containing Receptor 1 ELISA Kit (orb1213223)

Description	Human Epithelial Discoidin Domain-containing Receptor 1, DDR1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	DDR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	101,128
Uniprot ID	Q08345
Solubility (25°C)	24ng/ml
Sensitivity	0.028ng/ml
Entrez	780
Expiration Date	Please enquire.