

Product Datasheet

Human Ephrin Type-A Receptor 2 ELISA Kit (orb1213227)

Description	Human Ephrin Type-A Receptor 2, EPHA2 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	EPHA2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	108,266
Uniprot ID	P29317
Solubility (25°C)	24ng/ml
Sensitivity	0.026ng/ml
Entrez	1969
Expiration Date	Please enquire.