

Product Datasheet

Human Endoplasmic Reticulum Resident Protein 44 ELISA Kit (orb1213232)

Description	Human Endoplasmic Reticulum Resident Protein 44, ERP44 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	ERP44
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,971
Uniprot ID	Q9BS26
Solubility (25°C)	4800ng/L
Sensitivity	9.33ng/L
Entrez	23071
Expiration Date	Please enquire.