

Product Datasheet

Human Emilin-1 ELISA Kit (orb1213235)

Description	Human Emilin-1, EMILIN1 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	EMILIN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	106,695
Uniprot ID	Q9Y6C2
Solubility (25°C)	48ng/ml
Sensitivity	0.048ng/ml
Entrez	11117
Expiration Date	Please enquire.