

Product Datasheet

Human DNA-binding Protein Inhibitor Id-2 ELISA Kit (orb1213242)

Description	Human DNA-binding Protein Inhibitor Id-2, ID2 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	ID2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,917
Uniprot ID	Q02363
Solubility (25°C)	3200ng/L
Sensitivity	9.48ng/L
Entrez	3398
Expiration Date	Please enquire.