

Product Datasheet

Human DNA-binding Protein Inhibitor Id-1 ELISA Kit (orb1213243)

Description	Human DNA-binding Protein Inhibitor Id-1, ID1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	ID1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,133
Uniprot ID	P41134
Solubility (25°C)	4800ng/L
Sensitivity	10.11ng/L
Entrez	3397
Expiration Date	Please enquire.