

Product Datasheet

Human Disintegrin and Metalloproteinase Domain-containing Protein 28 ELISA Kit (orb1213251)

Description	Human Disintegrin and Metalloproteinase Domain-containing Protein 28, ADAM28 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	ADAM28
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	87,148
Uniprot ID	Q9UKQ2
Solubility (25°C)	24ng/ml
Sensitivity	0.018ng/ml
Entrez	10863
Expiration Date	Please enquire.