

Product Datasheet

Human Dickkopf-related Protein 2 ELISA Kit (orb1213256)

Description	Human Dickkopf-related Protein 2, DKK2 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	DKK2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,447
Uniprot ID	Q9UBU2
Solubility (25°C)	48ng/ml
Sensitivity	0.026ng/ml
Entrez	27123
Expiration Date	Please enquire.