

Product Datasheet

Human Deoxyribonuclease Gamma ELISA Kit (orb1213260)

Description	Human Deoxyribonuclease Gamma, DNASE1L3 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	DNASE1L3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,504
Uniprot ID	Q13609
Solubility (25°C)	3200ng/L
Sensitivity	7.55ng/L
Entrez	1776
Expiration Date	Please enquire.