

Product Datasheet

Human Death Domain-associated Protein 6 ELISA Kit (orb1213263)

Description	Human Death Domain-associated Protein 6, DAXX ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	DAXX
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	81,373
Uniprot ID	Q9UER7
Solubility (25°C)	6400ng/L
Sensitivity	9.38ng/L
Entrez	1616
Expiration Date	Please enquire.