

Product Datasheet

Human Cytokine Receptor-like Factor 1 ELISA Kit (orb1213266)

Description	Human Cytokine Receptor-like Factor 1, CRLF1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	CRLF1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,302
Uniprot ID	075462
Solubility (25°C)	24ng/ml
Sensitivity	0.019ng/ml
Entrez	9244
Expiration Date	Please enquire.