

Product Datasheet

Human Cytochrome P450 2C9 ELISA Kit (orb1213269)

Description	Human Cytochrome P450 2C9, CYP2C9 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CYP2C9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	55,628
Uniprot ID	P11712
Solubility (25°C)	3200ng/L
Sensitivity	7.59ng/L
Entrez	1559
Expiration Date	Please enquire.