

Product Datasheet

Human Cytochrome P450 26A1 ELISA Kit (orb1213270)

Description	Human Cytochrome P450 26A1, CYP26A1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	CYP26A1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	56,199
Uniprot ID	043174
Solubility (25°C)	2400ng/L
Sensitivity	4.9ng/L
Entrez	1592
Expiration Date	Please enquire.