

Product Datasheet

Human Cysteine Dioxygenase Type 1 ELISA Kit (orb1213273)

Description	Human Cysteine Dioxygenase Type 1, CDO1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	CDO1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,972
Uniprot ID	Q16878
Solubility (25°C)	4800ng/L
Sensitivity	10.04ng/L
Entrez	1036
Expiration Date	Please enquire.