

Product Datasheet

Human C-X-C Chemokine Receptor Type 6 ELISA Kit (orb1213276)

Description	Human C-X-C Chemokine Receptor Type 6, CXCR6 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	CXCR6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,280
Uniprot ID	000574
Solubility (25°C)	4800ng/L
Sensitivity	9.84ng/L
Entrez	10663
Expiration Date	Please enquire.